




Bringing Nature and Technology Together

Standard Reference Solution @ 25°C

*for calibration of
pH Meters and
electrode storage*

INSTRUCTIONS

Pour an adequate quantity of reference solution into a clean container. Immerse the pH electrode in pH reference solution for 1 minute before calibrating. To extend the life of your pH meter, keep electrode immersed in a small cup of pH reference solution when the meter is turned off. Rinse the electrode with distilled water after each pH reading to avoid contamination of the storage reference solution. Replace storage solution often. Never pour used solution back into this bottle.



**KEEP OUT OF REACH
OF CHILDREN**

Manufactured by
 **GH**
GENERAL HYDROPONICS
PO Box 1576, Sebastopol, CA 95473, USA
genhydro.com



946 ml (1 Quart)

s072506am